CICbiomaGUNE

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Vanessa Gómez-Vallejo · 943 005 330 · vgomez@cicbiomagune.es

CICbiomaGUNE

MEMBER OF
BASQUE RESEARCH
& TECHNOLOGY ALLIANCE

Paseo Miramón 194
Parque Tecnológico de Gipuzkoa
20014 · Donostia - San Sebastián
Gipuzkoa

Tel.: +34 943 005 300

www.cicbiomagune.es

CIC biomaGUNE

MEMBER OF
BASQUE RESEARCH
& TECHNOLOGY ALLIANCE

SERVICES CATALOG

2026

Acknowledgments

Centre accredited within the Basque Network of Science, Technology and Innovation by the Basque Government.

CIC biomaGUNE's R+D+i Management System is certified by SGS in terms of compliance with the ISO 56001 standard.

The CIC biomaGUNE Molecular and Functional Imaging Platform is part of the network of Singular Scientific and Technical Infrastructures (ICTS) of the Ministry of Science and Innovation.

The Animal Unit Platform is certified by the Association for the Evaluation and Accreditation of Laboratory Animal Care *AAALAC International*.

Centre accredited by the European Commission with the *HRS4R* Quality Seal for Human Resources Strategy for Research.

1. INTRODUCTION	4	
2. CONSULTING, TRAINING AND REPORTING	5	
3. SURFACE ANALYSIS AND COATING	6	
4. MASS SPECTROMETRY	9	
5. OPTICAL MICROSCOPY	11	
6. FLOW CYTOMETRY	14	
7. OPTICAL SPECTROSCOPY	15	
8. CALORIMETRY	17	
9. ELECTRON MICROSCOPY	18	
10. NUCLEAR MAGNETIC RESONANCE	20	
11. MOLECULAR AND FUNCTIONAL IMAGING UNIT		
11.1 ANIMAL FACILITY ANIMAL WELFARE	22	
11.2 MAGNETIC RESONANCE IMAGING	23	
11.3 RADIOIMAGING AND IMAGE ANALYSIS	24	
11.3 RADIOCHEMISTRY	26	

General Terms and Conditions of Sale

<https://www.cicbiomagune.es/general-terms-and-conditions-sale>

Posted on January 26, 2026

1. INTRODUCTION

One of the main elements of CIC biomaGUNE's mission is to improve the competitiveness of companies and support the development of the industrial fabric and scientific knowledge through the transfer of knowledge and the provision of our technical infrastructures.

This document contains a detailed description of the standardized services provided by the different technological platforms to support research in CIC biomaGUNE.

Contact the person responsible for the service you wish to contract.

All services will be carried out upon acceptance of a quote.

Possible variations in prices (upward or downward) on the requested services will be verified by the team of specialists depending on the type of analysis carried out in a new quote.

The prices indicated are presented in euros (€) and do not include value added tax (VAT). There are two general rates depending on the activity and nature of the applicant entity:

Rate A applies to non-profit R+D entities and entities associated with CIC biomaGUNE.

Rate B applies to all other organizations.

The prices indicated in this catalog are valid until January 31, 2027.

Visit the https://www.cicbiomagune.es/technology_offer website for current rates.

943 005 300 · external.services@cicbiomagune.es

2. CONSULTING, TRAINING & REPORTING

2.1 . SCIENTIFIC-TECHNICAL CONSULTANCY AND TRAINING

Scientific-technical consulting services are offered in different areas of specialization.

The technical team of the technological platforms offers advanced theoretical and practical training in scientific equipment in various areas of interest such as liquid chromatography, mass spectrometry, electron microscopy, optical microscopy, fluorescence microscopy, radiological protection, X-ray-induced photoelectron spectroscopy, etc.


2.2. PREPARATION OF REPORTS

The results obtained will be sent electronically to the address provided by the client. In general, the services offered include the cost of the personnel in charge of carrying out the service.

The preparation of reports derived from the analyses carried out, which entail the interpretation of the results, bibliographic search or comparative analyses between the samples, must be expressly requested by the client from the technical team of the CIC biomaGUNE technological platforms. The associated fee will be estimated in a budget specifically for the preparation of the report.

Contact us at the address below to request specific training according to your needs, or to evaluate other of our services.

943 005 300 · external.services@icbiomagune.es

SURFACE ANALYSIS AND COATING

3. SURFACE ANALYSIS AND COATING

3.1. DESCRIPTION OF SERVICES

Atomic force microscopy. It allows the topographic study of surfaces in samples dried in the air or immersed in liquid. From the topographic image with atomic resolution, the presence of impurities in materials can be determined. In addition, operating in force mode, indentation studies can be carried out, which allow quantitative maps to be obtained with mechanical properties such as the moduli of elasticity, deformation-dissipation or adhesion.

Applications: determination of surface roughness, grain size or size of nanoparticles on surfaces and substrates, nano-mechanical properties of surfaces. Applicable to the study of different substrates such as films, polymers, proteins, viruses, bacteria, biomolecules, etc.

Elemental composition analysis by X-ray-induced photoelectron spectroscopy. It is a surface analysis technique that allows the composition of the most superficial layer of a material to be determined, both qualitatively and semi-quantitatively. It also allows the chemical state of the elements present on the surface of the material to be determined. Sample analyses include a results report containing the general spectrum, the high-resolution spectra of the requested regions, the applied mathematical adjustment, and the calculation of the relative composition of the oxidation states for each region.

Applications: study of corrosion of materials, evaluation of substrate stability, determination of impurities, study of variations in the composition of materials as a function of the treatment time with argon plasma.

Physical vapor deposition. It allows the coating of surfaces (organic or inorganic) with mono or multiatomic layers with thicknesses ranging from one nanometer to several microns by means of magnetron-assisted *sputtering*.

Applications: metallization of samples to generate conductive surfaces, coating with semiconductors, bio-compatible coatings with titanium, improvement of the mechanical properties of materials, etc.

The following whites are available: Aluminum, Chromium, Copper, Hafnium, Niobium, Silicon, Tantalus, Titanium, Gold and Zinc. These elements can be deposited in pure form and in their corresponding carbide, nitride or pinworm derivatives. The feasibility of deposition may be evaluated with elements that do not appear in the list.


3.3. EQUIPMENT

Atomic force microscope, AFM – Veeco Bruker Multimode 8 HR-U. It is equipped with the **PeakForce QNM** module. It has a stand for measuring samples in liquid. Viscoelasticity

studies can be performed at frequencies between 0.1 Hz and 4 kHz. The maximum sample size is 1 cm². The system is controlled with the *Nanoscope Analysis 2.0 (Bruker) software*.

X-ray photoemission source scanning microscope, XPS – PHI VersaProbe III. It is equipped with a monochromatic X-ray source (K α aluminium at 1487 eV) that allows sweeping with a focused beam between 9 μ m and 200 μ m. The system can obtain secondary electron images and maps of the chemical state. 180° hemispherical analyzer and a detector with 128 channels. The motorized sample holder allows the center sample to be rotated symmetrically (ZalarTM rotation). Allows variable-angle XPS to be performed on samples up to 60 x 25 mm in size and temperatures between 500 and -120°C. It has an Argon ion source to neutralize the surface charge and perform depth profiles; it also has an Ultraviolet source to perform ultraviolet photoelectron spectrometry (UPS).

Magnetron-assisted spray vapor physical deposition coating system, PVD – AJA ATC 1800-UVH. It consists of an ultra-high vacuum chamber, and a power control system that with two radio frequency and direct current (RF/DC) power supplies, with the capacity to atomize and spray on the surface of any type of substrate. The rotary sample platform has a diameter of 9 cm and allows substrates to be heated to 600 °C.

3.4. FEE LIST AND CONTACT


Art. No.	Description · unit	A (€)	B (€)
SAF-AFM	AFM Bruker · h	93	107
SAF-AFM-TipA	Standard Tip for AFM · unit	52	60
SAF-AFM-TipB	Special Tip for AFM · unit	64	73
SAF-AFM-SA	AFM sample preparation · h	64	73
SAF-PVD-SP	PVD AJA ATC 1800 Sample Preparation	49	56
SAF-PVD-Au5	5 nm Gold · layer	19	22
SAF-PVD-Au100	100 nm Gold · layer	54	61
SAF-XPS-hour	XPS PHI Acquisition Measurement · h	109	124
SAF-XPS-Ins	XPS sample preparation and insertion time · h	64	73
SAF-Re	Report preparation · h	67	76

To: non-profit R+D entities, and entities associated with CIC biomaGUNE.

B: other organizations.

RESPONSIBLE FOR THE SERVICES OF SURFACE ANALYSIS AND COATING:

Desirè Di Silvio · Tel: 943 005 337 · ddisilvio@cicbiomagune.es

4. MASS SPECTROMETRY

4.1. DESCRIPTION OF SERVICES

Specialized technical advice for the development and planning of analysis of compounds or their derivatives.

Optimization of analyte determination by multiple characterization by combining multiparametric information from available equipment.

Determination of major components in liquid or solid matrices, with the possibility of making composition maps.

4.2. EQUIPMENT

Liquid chromatograph with high-resolution mass spectrometer, UHPLC/ESI-Q-ToF-MS – Waters Synapt XS. It is equipped with an *electrospray* ionization module and LockSpray interface to an orthogonal acceleration time-of-flight mass spectrometer (oa-TOF).

High Pressure Liquid Chromatograph, HPLC – Agilent LC 1260 Infinity II. It is equipped with the latest generation of DAD HS diode array detectors, a quaternary pump with an integrated degassing module and an autosampler with sample temperature control, programmable between 4 and 40 °C.

Inductively coupled plasma mass spectrometer, ICP-MS – Thermo Fisher iCAP-Q. It is equipped with a He Gas Kinetic Energy Discrimination Collision/Reaction Cell (KED). It is a robust analytical equipment, with high sensitivity and reproducibility that allows the quantification of most of the elements of the periodic table (metals and several non-metals), with detection limits in the ppb/ppt range.

Microwave digestion system, Berghof Speedwave Xpert. It allows microwave-assisted sample treatment to be carried out at high pressure and temperature in high-quality Teflon reactors, prior to ICP-MS analysis.

Matrix-assisted Laser Ionization/Desorption Tandem Mass Spectrometer, MALDI-TOF – neoflex. The Bruker neoflex is a benchtop MALDI-TOF/TOF mass spectrometer that allows the analysis of small molecules, peptides and intact proteins, with a mass range reaching approximately 300 kDa. The system offers high sensitivity for MS and MS/MS experiments and is specially designed for MALDI imaging of tissues and applications in top-down and bottom-up proteomics, lipidomics and spatial metabolomics.

MASS SPECTROMETRY

4.3. FEE LIST AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
MS-UPLC-MS	UPLC/ESI-Q-ToF-MS · h	127	146
MS-HPLC-DAD	HPLC Agilent LC 1260 Infinity II · h	64	74
MS-Sa	Mass Spectrometry Sample preparation · h	57	65
MS-MALDI-hour	MALDI-ToF Bruker · h	123	142
MS-MALDI-SR	MALDI-ToF Bruker · Standard Resolution · sample	19	22
MS-MALDI-HR	MALDI-ToF Bruker · High Resolution · sample	34	38
MS-ICP-Eq	ICP-MS iCAP-Q Calibration · service	90	103
MS-ICP-Sa	ICP-MS iCAP-Q Additional Samples · sample	13	14
MS-ACD-Opt	Digest optimization for ICP · sample	133	154
MS-AcD-6	Microwave Digestor (Up to 6 samples) · sample	72	83
MS-AcD-12	Microwave Digestor (Up to 12 samples) · sample	122	140
MS-Re	Report preparation · h	67	76

To: non-profit R+D entities, and entities associated with CIC biomaGUNE.

B: other organizations.

RESPONSIBLE FOR MASS SPECTROMETRY SERVICES:

Javier Calvo · 943 005 348 · jcalvo@cicbiomagune.es

5. OPTICAL MICROSCOPY

5.1. DESCRIPTION OF SERVICES

In vitro and ex vivo imaging of materials, biomaterials and biological systems, compositional analysis. It allows transmission, reflection and fluorescence images of biomaterials and biological samples with low or high resolution, both laterally and at depth.

Spectroscopic and chemical characterization of molecules and materials. Analysis and identification of organic, inorganic and biological molecules and materials, polymers, nanoparticles, etc., from solid, liquid and suspended samples.

Characterization of material composition, images and chemical profiles, three-dimensional chemical reconstruction of solid samples active in Raman with high spatial and spectral resolution.

5.2. EQUIPMENT

Fluorescence and multiphoton confocal microscope, Zeiss LSM 880 Airyscan and MaiTai DeepSee. It is equipped with five laser excitation sources from Ar, DPSS and HeNe operating at 405, 458, 488, 514, 561, 594 and 633 nm respectively. It has a multiphoton laser with variable excitation range from 690 to 1040 nm **MaiTai DeepSee**. It is equipped with two PMT photomultiplier detectors, one standard and the other with high sensitivity in the infrared range, as well as a GaAsP detector for detection in the visible range. The three spectral detectors can be used simultaneously in confocal imaging. In addition, it has an array of **Airyscan detectors** that allows measurements to be made with high lateral resolution of 150 nm and axial resolution of 400 nm, on channel 488. It also has a **BIG-2** detector module, consisting of two high-sensitivity GaAsP detectors for multiphoton measurements. Allows for fluorescence cross-correlation and correlation spectroscopy measurements.

Zeiss LSM 900 fluorescence confocal microscope. Fully automated, fully automated, inverted laser scanning fluorescence confocal microscope, Zeiss LSM 900. It enables simultaneous acquisition of high-resolution images in up to three fluorescence channels. It incorporates an AI Sample Finder system for rapid sample location and the execution of automated experimental routines. It allows the acquisition of Z stacks, time series and mosaics, as well as the performance of FRAP and FRET experiments.

It is equipped with excitation laser lines at 405, 488, 561 and 640 nm, two PMT detectors and a GaAsP detector for fluorescence, as well as an ESID detector for brightfield imaging. The system achieves a maximum resolution of 6144 × 6144 pixels and operates with fully licensed ZEN blue software, including advanced acquisition, co-location and 3D

OPTICAL MICROSCOPY

visualization tools. It has a wide range of EC Plan-Neofluar and Plan-Apochromat objectives for various applications in confocal microscopy.

Epifluorescence microscope, Zeiss Axio Cell Observer. It has a lighting system consisting of four Colibri LEDs (365, 470, 530 and 590 nm) and a Xe lamp. It is equipped with two AxioCam cameras: MRm3 and 305color. The microscope is equipped with an incubation booth for experiments with temperature, CO₂ and humidity control, which allows long-term *in vitro* experiments to be carried out. It also has a motorized plate for the characterization of tissues or samples in the millimeter range both by fluorescence and in brightfield (phase contrast, DIC and polarization). The system is controlled with Zen Blue software.

Raman confocal microscope, Renishaw inVia Reflex. It allows analysis of solids or substrates immersed in water, as well as liquid or suspended samples. The Leica microscope integrated into the system is equipped with objectives of different magnifications and a water immersion objective, a motorized scanning sample holder platform and a video camera. It is equipped with three laser excitation sources at 532, 633 and 785 nm (power 35, 17 and 180 mW). The spectrometer has two CCD cameras with diffraction gratings with 1800 and 1200 lines/mm that provide high spectral resolution (1 cm⁻¹) and another diffraction grating with 600 lines/mm optimized to measure at 532 and 633 nm with lower spectral resolution. In addition, it has polarization filters, StreamLine sweeping, sample accessory for macroscopic liquids. The system is controlled with Wire 4.4 software.

Raman confocal microscope, WITec alpha300 R. Allows the analysis of solids or substrates immersed in water. The integrated Zeiss microscope is equipped with objectives of different magnifications, a darkfield objective and two water-immersion objectives; a CMOS video camera and a topographic surface sensor that allows measurements on rough substrates while maintaining focus. It has four laser excitation sources coupled at 488, 532, 633 and 785 nm (power 55, 50, 30 and 75 mW), all equipped with TruePower, as well as two high-efficiency spectrometers. The visible-range spectrometer is equipped with 2 diffraction gratings for high and low spectral resolution (1800 and 300 lines/mm) and a CCD (EM) camera. The NIR spectrometer also has 2 diffraction gratings with 1200 and 300 lines/mm and a (BI) NIR CCD camera. The motorized scanning platform allows movement over an area of 50 x 50 mm² and is coupled to a piezoelectric scanner (100 x 100 x 20 μm) that allows confocal images to be obtained with a resolution close to the diffraction limit. The system is controlled with Control FIVE and Project FIVE software.

5.3. FEE LIST AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
OS-LSM900	Confocal Microscope Zeiss LSM 900 · h	106	121
OS-LSM880	Confocal Microscope Zeiss LSM 880 · h (Standard mode)	113	129
OS-LSM880-MP	Confocal System Zeiss LSM 880 · h (Multiphoton mode)	144	165
OS-CAO	Epifluorescence Microscope Zeiss Cell Axio Observer · h	102	117
OS-CRM-Ren	Renishaw inVia Reflex Confocal Raman Microscope · h	102	117
OS-CRM-B&W	Confocal Raman Microscope WITec Alpha 300R · h	102	117
OS-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomaGUNE.

B: other organizations.

Changes to the microscopes' preset settings may affect the price of the service.

THESE SERVICES ARE DISCOUNTED FOR FREQUENT USE:

- 10% from 20 hours per month
- 15% from 30 hours per month
- 20% from 40 hours per month

RESPONSIBLE FOR THE OPTICAL MICROSCOPY SERVICES:

Irantzu Llarena · 943 005 341 · illarena@cicbiomagune.es

Judith Langer · Tel: 943 005 341 · jlanger@cicbiomagune.es

FLOW CYTOMETRY

6. FLOW CYTOMETRY

6.1. DESCRIPTION OF SERVICES

Flow cytometry allows characterizing cellular systems and suspended materials in the micrometer range. It allows multiparametric analysis of cell populations or micromaterials according to the physicochemical characteristics of the sample, being able to evaluate 1000 events per second.

6.2. EQUIPMENT

Flow cytometer, FCM, MACSQuant Analyzer 10 Miltenyi Biotec. It has three excitation lasers: violet, blue and red (405, 488 and 640 nm). It has sample holders to measure on plates or tubes automatically. *MACSQuantify software* controls sample acquisition and enables data analysis.

6.3. FEE SCHEDULE AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
OS-FCM	Flow Cytometry · h	78	90
OS-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomaGUNE.

B: other organizations.

RESPONSIBLE FOR FLOW CYTOMETRY SERVICES:

Irantzu Llarena · 943 005 341 · illarena@cicbiomagune.es

7. OPTICAL SPECTROSCOPY

7.1. DESCRIPTION OF SERVICES

Measurement of the optical properties of compounds and biomolecules in solution, suspension or solid state. Measurement of the extinction or transmission of solutions, dispersions, or other partially transparent materials. Determination of wavelength and intensity of fluorescence emission, polarized fluorescence phosphorescence, emission lifetime, absorbance, transmission, circular and linear dichroism signal, etc.

7.2. EQUIPMENT

Spectrofluorometer, Edinburgh FS5. It allows the acquisition of emission, excitation and synchronous spectra from ultraviolet to near-infrared UV-NIR (between 200 and 1650 nm). It also allows for kinetic profile recording at multiple wavelengths, polarization measurements, and fluorescence anisotropy measurements on liquid, thin-film, and **solid samples**. It has an integrative sphere that allows measurements to be made to determine the quantum fluorescence yield and acquire diffuse reflectance spectra. The system allows the characterization of fluorescence and phosphorescence lifetimes, for this it uses high-intensity pulsed LEDs (340 and 405nm) and xenon pulsed lamp respectively. The system is controlled by Fluoracle *software*.

Circular dichroism spectropolarimeter, CD – Jasco J1500. It allows absorbance, circular and linear dichroism measurements from ultraviolet to near-infrared UV-NIR (180 and 1600 nm) and fluorescence measurements. It has a Xenon light source (power 150 W), Peltier temperature control and integrating sphere for measurements of samples with high light dispersion, solids, either by reflection or by diffuse reflectance. The system is coupled to an ATS-429L automated titration system. The system is controlled with Spectra *Manager II software*.

Fourier transform infrared spectrophotometer, FT-IR – Bruker Invenio-X. It allows solids light absorption analyses on solid, liquid or gel samples, in the transmission range between 30 and 15,500 cm^{-1} with a high spectral resolution and has an accessory with a simple reflection diamond monolith to measure the total attenuated reflectance (ATR). The equipment has a system for configuring beam *splitters* and detectors: CaF₂/InGaAs, KBr/DLaTGS and Si/DLaTGS/PE. The system is controlled by the *OPUS 8 software*. A hydraulic press is available for compression preparation of tablets.

Absorbance spectrophotometer, UV-Vis-nIR – Varian Cary 5000. It allows light absorption/extinction analyses to be performed on solid, liquid, suspended or thin film samples, in the range between 175 and 3300 nm. The measurement system is dual-beam and has a high spectral resolution (up to 0.01 nm in UV-Vis and 0.04 nm in nIR) thanks to its PMT (UV-Vis) and PbS (nIR) detectors. The spectral parameters can be adjusted

OPTICAL SPECTROSCOPY

individually for both ranges. It has a multi-cell stand that allows temperature-controlled measurements to be made within the range between 0 and 90 °C, via the *CaryWinUV* software.

Raman spectrometer, B&W Tek i-Raman 785. It allows analysis of macroscopic samples and materials in the form of liquids, gels, powders and solids. The spectrometer has a laser at 785 nm (power 250 mW) as the excitation source and operates in the spectral range from 175 to 3200 cm^{-1} with spectral resolution of 3.5 cm^{-1} . The light emitted and the reception of the signal are carried out by means of an optical fibre connected to the equipment, which has a CCD detector. Signal acquisition is controlled with *BWSpec4* software.

Dynamic light scattering, DLS Zetasizer Ultra Red Label. It allows measuring the molecular weight (from 1 kDa to 20 MDa), the hydrodynamic size of particles or molecules with a size in the range between 0.3 nm and 10 μm , determining the concentration (10^8 – 10^{12} particles mL^{-1}) and the surface charge of particles between 4 nm and 100 μm . The instrument is equipped with a 10 mW 633 nm laser, fluorescence filter, polarizer and photodiode detector. Signal acquisition and temperature control in the range of 4 to 120 °C is controlled with the *ZS Xplorer* software.

7.3. FEE LIST AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
OS-FS55	Fluorescence Spectrometer Edinburgh FS5 · h	67	76
OS-CDS	Circular Dichroism Spectrometer Jasco J-1500 · h	69	79
OS-FTIR-Sa	FT-IR Spectrometer Measurement · sample	23	26
OS-UV-nIR	UV-Vis-nIR Spectrophotometer Varian Cary 5000 · h	57	65
OS-RS	Raman Spectrometer B&W i-Raman 785 · h	57	65
OS-DLS	DLS Dinamic Light Scattering Zetasizer* · h	67	76
OS-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomaGUNE.

B: other organizations.

*The cost of consumables will be estimated according to the study.

RESPONSIBLE FOR THE OPTICAL SPECTROSCOPY SERVICES:

Irantzu Llarena · 943 005 341 · illarena@bicbiomagune.es

Judith Langer · Tel: 943 005 341 · jlanger@bicbiomagune.es

8. CALORIMETRY

8.1. DESCRIPTION OF SERVICES

Analysis and determination of enthalpies and thermal capacities of molecules and biomolecules in different solvents and mixtures by calorimetric methods: protein binding affinities, etc.

8.2. EQUIPMENT

Differential Scanning Calorimeter, DSC³⁺ – Mettler Toledo. It allows measuring the heat capacities of macromolecules. It is used to analyze and study materials such as thermoplastics, thermosets, elastomers, composite materials, metals and alloys, adhesives, food products, pharmaceuticals, and chemicals. The cooling/heating rate from 0.002 to 50/300 °K·min⁻¹ working in a temperature range of -50 to 500 °C.

Microscale thermophoresis, MST – Nanotemper Monolith. It allows the measurement of protein dissociation constants (*K_d*) (101-107 Da) in aqueous solution or cell lysates by adding serial dilutions of the ligand. It has two red and blue fluorescent channels. It allows the analysis of protein concentrations between 1 nM and 500 mM working in a temperature range between 20 and 40 °C with a sample volume of 4 µL.

8.3. FEE SCHEDULE AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
OS-DSC	Mettler DSC ³⁺ * Differential Scanning Calorimeter · h	19	22
OS-CM-Sa	Sample Preparation for Calorimeter Measurements · sample	33	38
OS-MST	Microscale Thermophoresis System* · h	90	103
OS-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomAGUNE.

B: other organizations.

* The cost of consumables will be estimated according to the study.

RESPONSIBLE FOR CALORIMETRY SERVICES:

Irantzu Llarena · 943 005 341 · illarena@cicbiomagune.es

ELECTRON MICROSCOPY

9. ELECTRON MICROSCOPY

9.1. DESCRIPTION OF SERVICES

The Electron Microscopy Platform has a Transmission Electron Microscope (TEM) optimized for ultra-high-resolution studies with low electron dose, which is equipped with: a CMOS camera, detector systems for the STEM scanning mode (BF and HAADF) and a high sensitivity X-ray scattered energy detector (EDS). It also has a Scanning Electron Microscope (SEM) that is equipped with an EDS system.

The platform is equipped with a dedicated sample preparation room that includes an immersion freezer, high-pressure freezer, automatic freezing replacement robot and an ultramicrotome with an attachment to process samples in cryogenic conditions.

The platform's expertise focuses on the preparation of (cryo)samples and corresponding (cryo)microscopy studies on nanoparticle compounds and polymeric and biological or soft materials, which can be prepared and transported to microscopes.

9.2. EQUIPMENT

Transmission electron microscope, TEM – JEOL JEM-2100F (UHR). It is an ultra-high-resolution device equipped with a FEG (*Field Emission Gun*) electron beam with variable acceleration between 80 and 200 kV, a 4-megapixel **TVIPS TemCam-F216** CMOS camera and two detectors for STEM (BF & HAADF) in combination with an EDS **Oxford UltimMax** detector and has the **GATAN 626 cryo-transfer system**.

Scanning electron microscope, FE-SEM – JEOL JSM IT800. It has secondary *Through-the-lens* (TTL) detectors and a main camera detector that provides excellent sensitivity to BSE. This set of lenses is reinforced by a solid-state detector system and a HAADF **GEN5** annular detector to perform SEM-STEM **DEBEN UK**, by an **ULTIM EXTREME SEM EDS system** controlled with the **AZTEC software** from the company **Oxford Instruments**, and also has a cryo-SEM system type **PP3010** from the company **QUORUM**.

Vitrification system, FEI Vitrobot Mark IV. It allows the cryogenics of samples in an automatic and reproducible way. It is used to carry out the process of vitrification of ultrafine suspended samples deposited on microscopy grids.

High pressure freezing system, Leica EM HPM 100. It enables the vitrification of biological and industrial samples with a typical thickness of 300 µm and a maximum thickness of 600 µm in a 3 mm diameter transport system.

Low-temperature cryosubstitution and embedding system, Leica EM AFS2. It allows cryosubstitution techniques with progressive temperature reduction (PLT), low-temperature inclusion and polymerization of resins.

Ultramicrotome, Leica EM UC7. It facilitates the preparation of ultra-thin and semi-thin slices of biological and industrial samples. This allows the analysis of sections, for example, by electron microscopy, optical or fluorescence. It has a **Leica EM FC7** cold working accessory chamber with which cuts can be prepared at temperatures between -15 °C and -185 °C.

9.3. FEE LIST AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
EM-TEM	TEM JEOL · h	153	175
EM-TEM-STEM	TEM · STEM Mode · h	173	199
EM-TEM-EDS	TEM · EDS Mode · h	173	199
EM-TEM-Cryo	TEM · Cryogenic samples · h	186	213
EM-TEM-Sa	TEM · Sample preparation (up to 4 samples)	19	22
EM-FESEM	FESEM JEOL JSM IT800 · h	97	111
EM-FESEM-EDS	FESEM · EDS mode · h	109	125
EM-FESEM-Cryo	SEM · Cryogenic samples · h	119	137
EM-FEI	Plunge Freezer System FEI Vitrobot Mark IV · h	92	105
EM-UC7	Leica UC7 Ultramicrotome · h	82	95
EM-UC7-Semi	Semi-thin sectioning (200–700 nm) · sample	61	69
EM-UC7-Ultra	Semi-thin sectioning (50–200 nm) · sample	74	85
EM-HPM	High Pressure Cryogenic Freezer Leica HPM100 · h	99	112
EM-AFS	Automatic Freeze Substitution Leica AFS2 · h	399	459
EM-LS	Lump-Sum Consumables · h	21	24
EM-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomagUNE.

B: other organizations.

RESPONSIBLE FOR ELECTRON MICROSCOPY SERVICES:

Marco Möller · 943 005 325 · mmoller@cicbiomagune.es

NUCLEAR MAGNETIC RESONANCE

10. NUCLEAR MAGNETIC RESONANCE

10.1. DESCRIPTION OF SERVICES

Nuclear Magnetic Resonance (NMR) is a very versatile spectroscopic technique widely used in many branches of chemistry, materials physics, biology or medicine. The most common applications include structural analysis, purity evaluation, structural elucidation, kinetic studies, identification and quantification of organic compounds, as well as metabolism studies in urine or blood plasma.

10.2. EQUIPMENT

Vertical nuclear magnetic resonance spectrometer, NMR – Bruker AVANCE III 500 MHz. The NMR equipment has a 5 mm BBI (inverse) Z-gradient automatic frequency tuning (ATM) probe for $^1\text{H}/^{19}\text{F}$ for standard sample analysis and allows measurements to be made at different temperatures. The platform has an HRMAS (*High Resolution Magic-Angle Spinning*) probe with three channels $^1\text{H}/^{13}\text{C}/^{31}\text{P}$ for the study of intact tissues or semi-solid samples.

Standard experiments:

- 1 dimension: ^1H , ^{13}C , ^{19}F , ^{29}Si , ^{31}P , ^{195}Pt , etc. Consult to measure other nuclides.
- 2 dimensions: COSY, TOCSY, HSQC, HMBC, NOESY, DOSY, etc.
- HRMAS (High Resolution Magnetic Resonance with Magic Angle).
- We adjust the experimental design to the need of your system.

Time-Domain Nuclear Magnetic Resonance Spectrometer, TD-NMR – Bruker Minispec mq60. It allows to determine the relaxation times, especially longitudinal and transverse relaxation times of MR (described by the T1 and T2 time constants, respectively) of compounds in solution or particles in suspension. Enables MRI contrast agent studies in a magnetic field of 1.4 T (60 MHz) and 37 °C

10.3. FEE SCHEDULE AND CONTACT

Art. No.	Description · unit	A (€)	B (€)
NMR-Sa	Sample preparation (Solvent & tube) · sample	12	14
NMR-Exp	NMR Experiment (First 15 min) · sample	37	42
NMR-hour	Additional NMR measurement time · h	2	2.3
NMR-Sa-MAS	Sample preparation (Solvent & probe) HRMAS · sample	30	35
NMR-SU-MAS	NMR Experiment HRMAS (First hour) · h	210	241
NMR-TD	NMR-TD Bruker Minispec mq60 · h	45	50
NME-Re	Report preparation · h	67	76

A: non-profit R+D entities and entities associated with CIC biomaGUNE.


B: other organizations.

THESE SERVICES ARE DISCOUNTED FOR FREQUENT USE:

- 10% from 40 measures per month
- 20% from 80 measures per month
- 30% from 120 measures per month

HEAD OF NUCLEAR MAGNETIC RESONANCE SERVICES:

Daniel Padró · 943 005 304 · dpadro@cicbiomagune.es

MOLECULAR AND FUNCTIONAL IMAGING UNIT

11. MOLECULAR AND FUNCTIONAL IMAGING UNIT

This 900 m² integrated facility offers state-of-the-art equipment for the synthesis (biomedical cyclotron + radiochemistry laboratory) and manipulation of compounds and nanosystems labeled with at least one radioisotope, for the realization of images by Positron Emission Computed Tomography (PET), Single Photon Emission Computed Tomography (SPECT) fused with X-ray Computed Tomography (CT) and Magnetic Resonance Imaging (MRI) images. in rodents housed in the Animal Facility Platform.

The unit is designed to perform *in vivo* studies from various imaging modalities on the same animal. The combination of these non-invasive imaging techniques facilitates multimodal approaches to biological, physiological and medical problems, obtaining images with functional and anatomical information. In addition, the unit is equipped with an autoradiography system for high-resolution nuclear imaging in biological tissues.

The Molecular and Functional Imaging Unit of CIC biomaGUNE is part of the Distributed Network of Biomedical Imaging (Singular Scientific and Technical Infrastructure, ICTS) formed by the nuclear imaging centers: BioImac (UCM), Imaging La Fe (UPV) and TRIMA (CNIC). Access through the network is via the www.redib.net website.

The services offered by each of the platforms that make up the Molecular and Functional Imaging Unit are described below.

11.1. ANIMAL FACILITY

The Animal Unit Platform is accredited by the AAALAC, has the capacity to house up to 1492 rodents (960 mice and 288 rats) and carries out comprehensive monitoring of animal welfare during procedures.

Comprehensive document management services are offered for licenses (drafting and application for animal experimentation projects specialized in molecular imaging); arrangements for the purchase, admission and shipment of animals for experimentation; as well as the maintenance, general care, cures and special diets of the animals.

We have qualified personnel in the administration of substances by routes: intravenous, intraperitoneal, subcutaneous, intracranial, intra-arterial or intramuscular, among others. Inoculation of type I tumor cells can be performed; general surgery and post-surgical follow-up; administration of analgesia and anesthesia, nanoparticles, radiotracers, MRI contrast agents, etc. We perform *in vivo* and *ex vivo* sampling, exsanguination, euthanasia, gross necropsy, including sending organs and formalin samples, etc.

The main functions of the Animal Unit Platform are:

- Ensure compliance with legal and ethical requirements on the use of animals for research and other scientific purposes.
- Provide care and monitor the health and welfare of laboratory animals of the species kept.
- Provide researchers with the advice and equipment needed to conduct their research on laboratory animals.

RESPONSIBLE FOR THE SERVICES OF THE ANIMAL UNIT PLATFORM:

Ainhoa Cano · 943 005 339 · acano@icbiomagune.es

11.2. MAGNETIC RESONANCE IMAGING

DESCRIPTION OF SERVICES

The Magnetic Resonance Imaging (MRI) Platform offers the instrumentation and expertise needed to perform a broad spectrum of magnetic resonance (MRI) imaging and spectroscopy studies applied to biological samples including small animals (rats, mice), tissue samples or extracts, and cell cultures. The facility has an operating room equipped for the preparation of animals for MRI exams; the implementation of surgical models and a wide variety of auxiliary equipment such as anaesthesia systems, physiological monitoring equipment compatible with MRI, infusion pumps, thermal regulation systems, etc.

The instruments are ideal for the following types of experiments:

- High-resolution 2D and/or 3D anatomical images.
- Diffusion and imaging of diffusion tensor and tractography.
- Flow imaging (angiography and tissue perfusion).
- Cardiac imaging.
- Chemical displacement imaging and localized spectroscopy.
- Development or implementation of new applications.

EQUIPMENT

11.7 T MRI Horizontal Spectrometer, MRI – Bruker BioSpec 117/16 USR. It has a cavity of 16 cm in internal diameter. It is equipped with a 9 cm (750 mT/m) gradient coil system and connected with a **Bruker AVANCE III** console (500 MHz) configured with 4 transmission channels and 8 reception channels.

7 T horizontal MRI spectrometer, MRI – Bruker BioSpec 70/30 USR. It has a cavity of 30 cm in internal diameter. It is equipped with two gradient coil systems: 20 cm (200 mT/m) and 12 cm (400 mT/m) and connected with a **Bruker AVANCE III** console (300 MHz) configured with 2 transmission channels and 4 reception channels.

MOLECULAR AND FUNCTIONAL IMAGING UNIT

Coils of different volume/surface sizes are available, optimized for brain, heart and body imaging for rats and mice, including *arrays* for the detection of other nuclei such as ^{13}C , ^{19}F and ^{31}P . The necessary equipment is also available to record the physiological state of the animal within the equipment, which allows the synchronization of the data acquisition of the scanner with the cardiac and respiratory movement.

FEE LIST AND CONTACT

Art. No.	Description	
MRI-7T	MRI Bruker Biospec Scanner 7 T (70/30 USR) · h	Request a Quotation
MRI-11T	MRI Bruker Biospec Scanner 11.7 T (117/16 USR) · h	Request a Quotation

RESPONSIBLE FOR THE SERVICES OF THE MAGNETIC RESONANCE IMAGING PLATFORM:

Daniel Padró · 943 005 304 · dpadro@cicbiomagune.es

11.3. RADIO IMAGING AND IMAGE ANALYSIS

DESCRIPTION OF SERVICES

The Radioimaging and Image Analysis Platform offers state-of-the-art preclinical imaging resources in Positron Emission Tomography (PET), Single Photon Emission Computed Tomography (SPECT) fused with X-ray Computed Tomography (CT) images for small animals. More specifically, it is equipped with a trimodal PET-SPECT-CT **MOLECUBES γ , β and X-CUBES** imaging system and a hybrid PET-CT **GE eXplore Vista-CT** system. The platform is designed to perform *in vivo* studies of several imaging modalities alternately on the same animal. The combination of these powerful non-invasive imaging techniques facilitates multimodal approaches to biological, physiological and medical problems, obtaining images with important functional and anatomical information. The platform is also equipped with an **Amersham Typhoon Biomolecular Imager Cytiva autoradiography system**, which allows high-resolution nuclear images of biological tissues to be obtained.

The Image Analysis Platform processes and analyzes all the multimodal images obtained within the Molecular and Functional Imaging Unit (PET, SPECT, CT and MRI). To do this, different workstations (PC, MAC) and *software* in various programming languages are available to develop automatic quantification and segmentation algorithms.

EQUIPMENT

Trimodal imaging system, PET, SPECT and CT - MOLECUBES γ , β and X-CUBES. These devices work separately allowing multimodal images to be obtained simultaneously. It is composed of three independent modules:

γ -CUBE. High-resolution, high-sensitivity SPECT equipment that allows the obtaining of full-body static images of rats and mice.

MOLECULAR AND FUNCTIONAL IMAGING UNIT

β -CUBE. High-resolution, high-sensitivity PET equipment with a wide axial field of view of 12 cm in a bed position, with which dynamic full-body acquisitions can be performed.

X-CUBE. High-performance CT equipment that allows high-resolution anatomical images (50 μ m) of rats and mice to be obtained in a short time.

Specific features of the system include: full-body mouse and rat imaging (PET, SPECT and CT); dynamic (PET) and static (PET and SPECT) studies; possibility of acquiring PET images synchronized with the heart rate of the animal *cardiac gated imaging*; possibility of acquiring CT images synchronized with the respiration of the animal *respiratory gated imaging*; possibility of acquiring PET images of up to three mice simultaneously; image reconstruction using the iterative 3DOSEM algorithm and scatter scatter correction and energy attenuation based on CT image.

Quadrimodal imaging system, PET, SPECT, CT and RO – MILabs VECTor. These devices work separately allowing multimodal images to be obtained simultaneously. It is composed of three independent modules:

Integrated multimodal system that combines PET, SPECT, CT and optical imaging techniques in a single, compact platform for small animal studies. It allows for dynamic and static full-body imaging of mice and rats.

It includes a high-resolution PET module for dynamic and multi-isotope studies; a SPECT system with resolution up to 0.5 mm and a wide energy range; a high-resolution CT (<25 μ m) with rapid full-body acquisitions; and an optical imaging module for fluorescence, bioluminescence, and Cherenkov effect with CT-based anatomical correction.

The system allows acquisitions synchronized with heart and respiratory rate, multi-isotope studies, simultaneous imaging of several animals, advanced reconstruction using iterative algorithms and export of images in standard DICOM and NIfTI formats.

Positron Emission Tomography Imaging System, PET-CT - GE eXplore Vista-CT. It provides high-sensitivity functional imaging along with anatomical imaging in a single instrument. The 7 cm cavity is suitable for housing mice and rats up to 400 g. With a total axial field of view of 15 cm (4.8 cm per position), this device allows both dynamic and static studies to be carried out. Specific features of the system include automatic bed movement during acquisition, FBP, 2DOSEM and 3DOSEM reconstruction methods, dispersion correction and CT image-based energy attenuation.

Autoradiography Equipment, Amersham Typhoon Biomolecular Imager Cytiva. It allows data to be acquired simultaneously from samples labeled with multiple radioisotopes, accurately discriminating the emissions of each isotope based on emission energy or

MOLECULAR AND FUNCTIONAL IMAGING UNIT

radioactive half-life. In a few hours, data can be acquired with the same quality and resolution as with traditional autoradiography films or phosphor screens. Decays are detected in real time, and the data visualization is continuously updated during acquisition, thus allowing optimal adjustment of exposure time.

FEE LIST AND CONTACT

Art. No.	Description	A (€)	B (€)
RIA-Mol	PET-SPECT-CT β γ , β and X-CUBES · h	Request a	Quotation
RIA-PET-GE	PET-CT GE eXplore Vista-CT · h	Request a	Quotation
RIA-ANI	Analysis of Nuclear Images · h	Request a	Quotation
RIA-BIA	Biomolecular Autoradiography Cytiva · h	Request a	Quotation
RIA-Re	Report preparation · h	Request a	Quotation

RESPONSIBLE FOR THE SERVICES OF THE RADIO IMAGE AND IMAGE ANALYSIS PLATFORM:

Unai Cossío · Tel: 943 005 316 · ucossio@cicbiomagune.es

11.4. RADIOCHEMISTRY

DESCRIPTION OF SERVICES

The Radiochemistry Platform offers the manufacture of radiotracers (radio-synthesis, purification, formulation and quality control) by incorporating at least one radioisotope into molecules or particles of interest.

The platform offers the production of the isotopes ^{11}C , ^{13}N , ^{15}O , ^{18}F , ^{64}Cu , ^{68}Ga and ^{89}Zr for use in the radiochemistry laboratory of CIC biomaGUNE, being able to label molecules of interest and be sent to other centers in the national territory.

The radiochemistry laboratory is equipped with a **Cyclone IBA 9/18** cyclotron with 7 targets (one of them for solid targets) and five leaded cells that house automatic and versatile synthesis modules, suitable to produce PET and SPECT radiotracers. The installation has modules specially designed to:

- Synthesis, purification, formulation and quality control of molecules with radioisotopes.
- Production and handling of radioisotopes: ^{11}C , ^{13}N , ^{15}O , ^{18}F , ^{64}Cu , ^{68}Ga and ^{89}Zr .
- Synthesis of ^{11}C CH₃I/ ^{11}C CH₃OTf from ^{11}C CO₂/ ^{11}C CH₄ and subsequent methylation reaction.
- ^{18}F -Fluorination of molecules of interest by nucleophilic or electrophilic substitution.
- Chelation reactions with radiometals: ^{68}Ga , ^{67}Ga , ^{64}Cu , ^{89}Zr , ^{177}Lu , etc.

In addition, the laboratory has four leaded cells suitable for the manual handling of radioactive substances, ideal for the preparation of single-dose or synthesis with small amounts of radioactivity.

MOLECULAR AND FUNCTIONAL IMAGING UNIT

Within the laboratory itself is the quality control area, which has state-of-the-art equipment to carry out the complete quality control of synthesized radiotracers, including the following equipment: radio-HPLC, radio-GC, radio-TLC and high-resolution multichannel detectors.

CATALOGUE OF RADIOTRACERS AND RADIONUCLIDES.

[¹¹ C]		[¹³ N]	
Choline	Amino acid · Cancer diagnosis	NH ₃	Ammonium · Perfusion and Blood Flow
DASB	Serotonin transporter	[¹⁸ F]	
FLB-457	Dopaminergic Receptor D2/D3	FAP	Fibroblast activation protein
Flumazenil	GABAA Receptors	FDG	F-Deoxyglucose · Glucose metabolism
Harmine	Enzyme Marker (MAO-A)	FLT	Tumor proliferation · Cancer diagnosis
Methionine	Amino acid · Cancer diagnosis	FMISO	Hypoxia
NNC-112	Dopaminergic Receptor D1	NaF	Fluoride · Bone alteration
GDP	Amyloid beta plaque	PE2I	Dopamine transporter
PK-11195	Benzodiazepine peripheral receptor (RPBs)	[⁶⁸ Ga]	
QNB	Cholinergic muscarinic receptor	DOTA	Somastostatin receptor
Raclopride	Dopaminergic Receptor D2	NODA	Somastostatin receptor

Non-encapsulated radionuclides approved at the Radioactive Facility: ³H, ¹¹C, ¹⁴C, ¹³N, ¹⁵O, ¹⁸F, ⁴⁵Ti, ⁵⁷Co, ⁶⁴Cu, ⁶⁷Ga, ⁶⁸Ga, ⁸⁹Zr, ⁸²Rb, ^{99m}Tc, ¹¹¹In, ¹²³I, ¹²⁴I, ¹²⁵I, ¹³¹I, ¹⁷⁷Lu, ²⁰¹Tl and ²²⁵Ac.

RESPONSIBLE FOR THE SERVICES OF THE RADIOCHEMISTRY PLATFORM:

Vanessa Gómez-Vallejo · 943 005 330 · vgomez@cicbiomagune.es

CICbiomaGUNE

MEMBER OF
BASQUE RESEARCH
& TECHNOLOGY ALLIANCE

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